


INSTRUMENT LAYOUT AND PERFORMANCE


First Proof of Concept Instrument at TOFWERK



N2ICP-001 Prove of Concept (POC) Instrument and Schematics

- MICAP plasma source (RADOM Corp.)
- Grounded sampler and skimmer
- Extraction lens
- 63 mm gate valve
- Ion Mirror
- Collision-/Reaction Cell
- Notch Filter
- TOFMS with R3000


Performance Characterization using Micro Droplet Generator (MDG) technology



- 100 ppb multi-element standards (Merck IV, REE's)
- Droplet volume: 142 pL
- Analyte per droplet: 14.2 fg
- 3 ms integration time
- Plasma on
 - Air (not dried)
 - N2 (Generator using PSA principle)

Video capture of Single droplet
Mass Spectrum and time dependent signals of single droplets


Passive Sampling by Venturi Effect



Air sampling using standard nebulizer

- Air sampling without pump in the sample flow path to the ion source (no GED required)
- Membrane pump to flush sample line
- Use of conventional Nebulizer
~0.15 L/min air sampling

Single particle signals, LOD's from filtered air



- Background: HEPA Filter
- Continuous measurement with 1 s integration time
- Background LOD determined from filtered aerosol
 - Log-Normal background fit accepting 5% false-positive/-negative

Signal from filtered and un-filtered outdoor air (left)
LOD calculation from filtered air signal (right)

DIRECT AIR SAMPLING – CASE STUDY

Are there Metals in the Air?



- Indoor Air
- TOFWERK shares workspace with solar industry
- Are there metals in the air?
- Are metals indicative for manufacturing activities?


TOFWERK headquarters in Thun (left)
Workspace shared with solar industry (right)

DIRECT AIR SAMPLING – PARTICLE INFORMATION

Application


2 days experiment

Signal with 1 s integration time



- Filtered and unfiltered Indoor Air
- Signal in [counts] converted to [fg] using droplet calibration.
- Sampled Analyte Mass per Time [fg/s]
- Large number of metals found in form of particulate matter
- Time traces indicate for different sources


Element ratios



- 1s integration time
- Ratios of mass per particle: Cr/Fe
- Stainless steel particles
high Cr content – soldering wire

Element ratio depending on particle mass

Average Concentration in ng/m³



Signal summed to 10 min intervals

- Sampled air volume per time [m³/s]
- Sampled Analyte Mass per Time [ng/s]
- Analyte concentration in sampled air [ng/m³]

Lower LOD than from single Particle Data


MAK values for Pb in air

0.1 mg/m³ - Pb and Pb-compounds (except Pb-Me₄ and Pb-Et₄)

0.05 mg/m³ - Pb-Me₄ and Pb-Et₄

Source: SUVA

Signal summed to 10 min integration time



- Signal measured at 1 s integration
- Smoothing of the Signal by summation to 10 min intervals
- Ability to follow dynamic events
- Washout after intense events?

Lower LOD than from single Particle Data

DIRECT AIR SAMPLING - SUMMARY

Application

MICAP-TOFMS provides

- Direct air sampling
- Single particle information [fg/particle]
- Detection of concentrations in [ng/m^3] range.
- Measurement intervals:
 - Milliseconds, minutes, days
- Ability to follow dynamic concentrations
- Indication for particle sources
 - Dynamic time dependent signal
 - Composition of individual particles





THANK YOU!!

Martin Tanner

m.tanner@tofwerk.com

27th ETH-Nanoparticles Conference (NPC-24)