

# New Technique for Online Measurements of Trace Metals in Ambient Particulate Matter (PM)

**Dongbin Wang<sup>1</sup>, Arian Saffari<sup>1</sup>, Martin M Shafer<sup>2</sup>, James J Schauer<sup>2</sup> and Constantinos Sioutas<sup>1\*</sup>**

*<sup>1</sup>Department of Civil and Environmental Engineering, University of Southern California, Los Angeles, California, USA*

*<sup>2</sup>Environmental Chemistry and Technology Program, University of Wisconsin-Madison, Madison, Wisconsin, USA*

19th ETH Conference on Combustion Generated Nanoparticles

# Motivation, Rationale and Objectives

- Transition elements and metals are important PM species found usually **in trace concentrations**
- Some are **highly toxic** (Fe, Cu, V, Cr, Mn, Ni etc)
- They are also important species for **PM source apportionment**
- Traditional PM sampling methodologies usually involve collection of particles onto substrates by **filtration or inertial impaction**
- Collected particles are **extracted using different types of solvents** and analyzed by various **analytical techniques** (ICP-MS, etc)
- **Drawbacks:**
  - **Long turn-around time** for processed results
  - **Long sampling intervals** for collection of sufficient mass for subsequent chemical analysis
  - **Changes in chemical speciation** (e.g. oxidation state shifts in redox active metals)
  - **Expensive analysis** (especially ICP-MS)

- **Advanced on-line aerosol sampling technologies** were developed in recent years to overcome these limitations (PILS etc)
- Most of these technologies are **inadequate for the measurement of trace level species** due to their **low signal-to-noise ratio**.

## Aerosol-Into-Liquid Collector:

**Aerosol-Into-Liquid Collector** (*Wang et al., AS&T, 2013*) provides highly concentrated slurry samples of ambient PM

### Advantages:

- High sampling flow rate (200 lpm)
- Provides **highly concentrated PM<sub>2.5</sub> and PM<sub>2.5-10</sub> slurry** samples
- Ability of **continuous stable operation** (at least 5 consecutive days)
- We will report on line measurement results for Cu, Fe, Cr and Mn
- Convert air pollution measurements into water pollution measurements




Figure 2 System schematic for collection efficiency tests

# Design of the New Impactor Sampler- Sampling flow 200 lpm (Wang et al AS&T 2013)




Figure 1 Schematic of impactor


# TABLE 1 Collection efficiency For Polystyrene Latex (PSL) particles and Black Carbon (BC)

Particle Size	Collection efficiency	Inlet loss	Wall loss	Fraction collected in suspension
PSL 42.5 nm	96.4%	4.5%	9.4%	82.4%
PSL 100 nm	90.3%	2.1%	7.5%	80.7%
PSL 300 nm	98.9%	1.1%	6.2%	91.4%
PSL 750 nm	95.1%	2.8%	8.1%	84.1%
PSL 1 $\mu\text{m}$	99.1%	1.7%	9.3%	88.0%
Ambient BC	94.3%	—	—	—

# Field evaluation

- Continuity tests/ Unattended Sampling

Ambient PM at USC were continuously collected over ~72 hrs period



- Number based collection efficiency = **98.8%**
- Gaps in OPS data due to software malfunction.

# On line sampler for near- continuous measurement of PM<sub>2.5</sub> and PM<sub>2.5-10</sub> Copper (Cu) Ion Selective Electrodes (ISE)

- Ion Selective Electrodes (ISEs) are routinely used for determining the concentration (*technically the activity*) of **selected ions in aqueous solutions**.
- Relatively **inexpensive, simple to use**, and mostly function linearly over many orders of magnitude of concentration.
- **Accurate measurement** (within a measurement error of 2% typically) in a **very short response time** (several seconds to several minutes).





## On Line PM-Bound Copper (Cu) Monitor






Figure 2 Schematic of **Cu on line measurement system**

# Comparison between Cupric ISE and ICP-MS analysis of off line filter samples

## Total Cu




## Water-soluble Cu



Overall very good agreement between cupric ISE and ICP-MS measurements were observed ( $R^2 > 0.87$  in both comparisons).


# Field evaluation

- Comparison between filter and slurry samples




The median slurry vs filter concentration ratio : 1.02 ± 0.12

- Continuous system operation




- Samples were collected/measured every 3-4 hours.
- Automated sample collection, measurement and system rinsing.
- Excellent performance for 5 consecutive days of unattended sampling.




# System configuration for on line measurements of Fe, Cr and Mn



## Sample Analysis Module




Species	Ligand
Cr (VI)	Diphenylcarbazide
Fe II (III)	Ferrozine
Mn II	Formaloxine (FAD)



Very good agreement between off-line ICPMS and online data

# Particle –Bound Reactive Oxygen Species (ROS) Measured by the On Line Sampler and off-line Filter and impactor samplers



Higher ROS with impactor/ on line sampler

Virtually identical ROS between samplers when insoluble species are removed

## PUBLICATIONS

1. Wang D., Kam W., Cheung K. and Sioutas C.\* “Development of a Two-Stage Virtual Impactor System for High Concentration Enrichment of Ultrafine, PM<sub>2.5</sub> and Coarse PM”. *Aerosol Science and Technology*, 47(3): 231-238, 2013
2. Wang D., Pakbin P., Shafer M.M., Schauer J.J. and Sioutas C.\* “Macrophage reactive oxygen species activity of water-soluble and insoluble fractions of ambient coarse, PM<sub>2.5</sub> and ultrafine particulate matter (PM) in Los Angeles “. *Atmospheric Environment*, 77:301-310, 2013
3. Wang D., Pakbin P., Saffari A., Shafer M.M, Schauer J.J and Sioutas C.\* “Development and Evaluation of a High-Volume Aerosol-Into-Liquid Collector for Fine and Ultrafine Particulate Matter”. *Aerosol Science and Technology*, 47(11):1226-1238, 2013
4. Wang D., Shafer M.M, Schauer J.J and Sioutas C.\* “Development of a technology for online measurement of water-soluble and total copper (Cu) in PM<sub>2.5</sub>.” *Aerosol Science and Technology*, 48(8), 864-874, 2014
5. Wang D., Shafer M.M, Schauer J.J and Sioutas C.\* “A new technique for online measurement of total and water-soluble copper (Cu) in coarse particulate matter (PM)” *Environmental Pollution*, 199, 227-234, 2015



# Acknowledgements

This study has been supported by :

- The **South Coast Air Quality Management District (AQMD)** through award number #11527
- **NIH grant number R01AI065617-13** to the University of Southern California (USC).
- The research described herein has not been subjected to the agency's required peer and policy review and therefore does not necessarily reflect the views of the agency, and no official endorsement should be inferred. Mention of trade names or commercial products does not constitute an endorsement or recommendation for use.