

A 2λ Photo-Thermal Aerosol Absorption Monitor: design, calibration and measurements of absorption enhancement due to coating

ETH nano Conference
Wed, 22 Jun 2022

L. Drinovec, U. Jagodič, L. Pirker, M. Škarabot, M. Kurtjak, K. Vidović, L. Ferrero, B. Visser, J. Röhrbein, E. Weingartner, D. M. Kalbermatter, K. Vasiliatou, T. Bühlmann, C. Pascale, T. Müller, A. Wiedensohler, **Griša Močnik**

[grisa.mocnik@ung.si](mailto:gresa.mocnik@ung.si)

Filter photometer biases








Photo-thermal interferometry

Particle absorb light




Pump beam intensity:
 $I=1 \text{ W/mm}^2$
Particle temperature:
 $\Delta T=3.3 \text{ K}$

Heat transferred to gas



Particle number:
 2500 cm^{-3}
Air temperature:
 $\Delta T=42 \mu\text{K}$
Refractive index:
 $\Delta n=47*10^{-12}$

Measurement of
interferometer
phase change




Light path length:
 $\Delta s=47*10^{-15} \text{ m}$
Interferometer phase:
 $\Delta \phi=0.15 \mu\text{Rad}$


Phase change is proportional to absorption!!

Instrument: PTAAM- 2λ

- Moosmüller, H. & Arnott, W. (1996). *Opt. Lett.*, 21, 438-440.
- Sedlacek, A.J. (2006). *Rev. Sci. Instrum.*, 77, 064903, 1-8.
- Visser et al. (2020). *Atmos. Meas. Tech.* 13, 7097–7111.
- Drinovec et al. (2022), *Atmos. Meas. Tech. Discuss.*, amt-2021-21

- Photo-Thermal Aerosol Absorption Monitor
- Mach-Zender interferometer, similar to Moosmuller, Arnott, Sedlacek and Visser
- Pump beam focused by **axicon** (patent EP 3492905)
- Simultaneous measurements at **532 and 1064 nm**





pump:532 nm & 1064 nm

Article



Atmospheric Measurement Techniques

<https://doi.org/10.5194/amt-2022-21>

<https://doi.org/10.5194/amt-2022-21>
© Author(s) 2022. This work is distributed under
the Creative Commons Attribution 4.0 License.



Abstract Discussion Metrics

22 Feb 2022

Review status: this preprint is currently under review for the journal AMT.

A dual-wavelength photothermal aerosol absorption monitor: design, calibration and performance

Luka Drinovec^{1,2,3}, Uroš Jagodič^{1,2}, Luka Pirker^{2,4}, Miha Škarabot², Mario Kurtjak⁵, Kristijan Vidović⁶, Luca Ferrero⁶, Bradley Visser⁶, Jannis Röhrbein⁸, Ernest Weingartner⁶, Daniel M. Kalbermatter⁹, Konstantina Vasilatou⁹, Tobias Bühlmann⁹, Celine Pascale⁹, Thomas Müller¹⁰, Alfred Wiedensohler¹⁰, and Griša Močnik^{1,2,3}

¹Haze Instruments d.o.o., Ljubljana, Slovenia

²Department of Condensed Matter Physics, Jozef Stefan Institute, Ljubljana, Slovenia

³Center for Atmospheric Research, University of Nova Gorica, Nova Gorica, Slovenia

⁴Faculty for mathematics and physics, University of Ljubljana, Ljubljana, Slovenia

⁵Advanced Materials Department, Jozef Stefan Institute, Ljubljana, Slovenia

⁶Department for Analytical Chemistry, National Institute of Chemistry, Ljubljana, Slovenia

⁷GEMMA center, University of Milano-Bicocca, Milano, Italy

⁸Institute for Sensors and Electronics, University of Applied Sciences Northwestern Switzerland, Windisch, Switzerland


⁹Federal Institute of Metrology METAS, Bern, Switzerland

¹⁰Leibniz Institute for Tropospheric Research, Leipzig, Germany

Received: 14 Jan 2022 – Accepted for review: 20 Feb 2022 – Discussion started: 22 Feb 2022

Photo-Thermal Aerosol Absorption Monitor




PTAAM 2 λ



The Eurostars Programme is powered by
EUREKA and the European Community



Photo-thermal signal



- Signal shape defined by the pump geometry
- Signal amplitude proportional to the absorption coefficient – lock-in amplifier

532 nm channel calibration

Calibration gas cylinder

Instrument response depends on:

- Pump beam intensity
- Overlap between pump and probe beams



-> instrument must be **calibrated**.



532 nm channel calibration:





- Using absorbing gas: **NO₂**
- Absorption cross-section = $1.47 \times 10^{-19} \text{ cm}^2$
- High uncertainty of NO₂ concentration < sub ppm range

Permeation NO₂ generator



1064 nm channel calibration

- lack of absorbing gasses in near-IR region
- calibration with **nigrosin** particles
 - absorbtion in VIS and IR
 - water soluble pigment -> forms nice **spherical** particles





Mie theory




Absorption coefficient

$$\text{babs}(1064 \text{ nm})/\text{babs}(532 \text{ nm}) \approx 0.075$$

Validation – AeroTox and MORE-RAIN campaign



Size distribution of aerosolized nigrosin




5-6% difference


in absorption coeff.
between
measurement and
Mie model

Results – linearity of response

NO₂



Nigrosin




Measurement uncertainty

Sources of uncertainty		Uncertainty	Components
A	NO ₂ amount fraction	2%	
B	Nigrosin refractive index	2%	
C	Mie calculation of $b_{\text{abs},1064 \text{ nm}}/b_{\text{abs},532 \text{ nm}}$	4%	
D	Scattering & absorbing gases	1%	
E	Stability of instrument	3%	
Combined uncertainties			
	$b_{\text{abs},532\text{nm}}$	4%	A, D, E
	$b_{\text{abs},1064\text{nm}}$	6%	A, B, C, D, E
	AAE	9%	B, C, D, E, ln

Low measurement uncertainty -> **reference method** for aerosol absorption measurements.

Lab soot – calibration of filter photometers

Campaign	Sample	Volume mode (nm)	Angstrom exp PTAAM (532 nm/1064 nm)	Angstrom exp AE33 (470 nm/950 nm)
Ljubljana	Diesel soot	130	1.05	1.31
AEROTOX	CAST soot	160	0.87	1.20
Ljubljana	Propane soot	400	0.86	1.03





$$C = \frac{b_{ATN(AE33)}}{b_{abs(PTAAM)}}$$

as expected from Mie theory for particles of this size





- C depends on the particle size
- higher C at 532 nm compared to 1064 nm
→ higher apparent Angstrom exp. in AE33

Filter photometers – SSA dependence



Absorption enhancement by coating

BC particles coated with α -pinene





Discrete
Dipole
Approximation

Liu et al., 2016

Kalbermater et al., 2022

Models are not adequate for description of absorption enhancement!

Ambient calibration of filter photometers




Conclusions

- **new photo-thermal aerosol absorption monitor**
- **calibrated to primary standards**
- **absorption coefficient at 532 nm & 1064 nm**
- **linear response & no artefacts**
- **low measurement uncertainty -> reference instrument for absorption**



**Thank you!
Questions?**

grisa.mocnik@ung.si



Swiss National Science Foundation (grant no. 200021_172649), EUROTARS programme (IMALA, grant no. 11386), Slovenian Research Agency grants P1-0385, P1-0099 and I-0033, the European Metrology Programme for Innovation and Research (EMPIR Black Carbon, EMPIR AeroTox), Paul Scherrer Institute, Teledyne API, EAS Envimet Analytical Systems GmbH, AMES d.o.o.

A note on the terminology

- **Black Carbon mass concentration** $BC \quad \mu\text{gm}^{-3}$
 - Optical meas. **equivalent** to mass **eBC** (filter photometers)
 - ...
- **Aerosol absorption coefficient** $b_{abs} \quad \text{Mm}^{-1}$
 - Optical measurements, in-situ, (in)direct

$$b_{abs} = BC \cdot MAC$$


C...c...calibration to determine $C \leftrightarrow MAC!$

Field: comparison w/ other filter photometers


Aerosol absorption coefficient b_{abs} Mm⁻¹

- need reference sample
- need **reference method to measure absorption**
- **Reference methods:**
 - extinction - scattering: low SSA
 - photoacoustics: non-coated particles, no VOCs or water
 - **photothermal interferometry:** sensitivity


Mach-Zehnder Interferometer



Mach-Zehnder Photothermal Interferometer




Quadrature Point



Highest **sensitivity** in φ change – the interferometer in **quadrature** ($I_1 = I_2 = \frac{I_0}{2}$), e.g. at $\varphi_o = \frac{\pi}{2}$.
 $\rightarrow \varphi$ change $\Delta\varphi = \varphi - \varphi_o = \frac{\Delta I}{I_o} = \frac{I_1 - I_2}{I_1 + I_2}$

Real part of the refractive index



- Real part for the UV-VIS region is taken from Bluvstein et al. (2017)
- Near infrared value was determined by Brewster angle measurement:
 $n(1064 \text{ nm}) = 1.848 \pm 0.005$



Nigrosin refractive index determination

Comparison of aquaous solution of nigrosin and thin film:

- solvent effect: +/-25%




$$n(532 \text{ nm}) = 1.62 + 0.223i$$

$$n(1064 \text{ nm}) = 1.73 + 0.0419i$$

Mie calculation

- Nigrosin aerosolisation using Topas ATM 226
- Particle size spectra measurements with SMPS
- Calculation of absorption using Mie routines in Matlab



MAC is stable for particles generated with 0.1 g/l nigrosin solution

$$\text{babs}(1064 \text{ nm})/\text{babs}(532 \text{ nm}) \approx 0.075$$